

Pre- Organized Hofmann Sets



Brief History



R a o u l H o f f m a n n












In 1938, Raoul Hoffmann, a German born surgeon living and working in Geneva, Switzerland, set himself the goal of designing a new way to reduce and fix broken bones. His objective was to devise a system that would utilize an alternative form of fracture treatment, without reliance on open surgery. The result was history's first fully-functional, mainstream External Fixation System

Pre-Organized Sets




By indications and as a set














Combinational external fixator(A type)

Combination	Name	NO.	Size	Picture
Ulna&Radiu spanning external fixator GX201001	Rod to rod coupling	201002-0505	φ5/5	
	4-hole pin clamp	201004-0534	φ5/3-4	
	Post	201009-0530	φ5/30°	
	Connecting rod	AL05-200	φ5×200	
	Bone screw	SS20301240	φ3×120	
Forearm GX201002	Rod to rod coupling	201002-0505	φ5/5	
	4-hole pin clamp	201004-0534	φ5/3-4	
	Post	201009-0500	φ5/0°	
	Connecting rod	AL05-200	φ5×200	
	Bone screw	SS20301240	φ3×120	
Distal radius GX201003	Rod to rod coupling	201002-0505	φ5/5	
	Pin to rod coupling	201001-0534	φ5/3-4	
	Peri-articular pin clamp	201003-0534	φ5/3-4	
	Connecting rod	CF05-200	φ5×200	
	Bone screw	SS20301240	φ3×120	
Wrist joint fixator GX201004	Wrist joint movement device	201010-05	φ5	
	Pin to rod coupling	201001-0534	φ5/3-4	
	Bone screw	SS20301240	φ3×120	
Elbow joint fixator GX201005	Elbow joint movement device	201011-05	φ5	
	Pin to rod coupling	201001-0534	φ5/3-4	
	Bone screw	SS20301240	φ3×120	
Mini-fixator (Finger&Plam) GX201006	Pin to rod coupling	201005-0302	φ3/2	
	Connecting rod	CF03-080	φ3×80	
	Bone screw	SS10200716	φ2×70	
Humeral shaft external fixator GX201007	Rod to rod coupling	201002-0505	φ5/5	
	4-hole pin clamp	201004-0534	φ5/3-4	
	Post	201009-0530	φ5/30°	
	Connecting rod	CF05-200	φ5×200	
	Bone screw	SS20401240	φ4×120	
Proximal humerus fixator GX201008	Rod to rod coupling	201002-0505	φ5/5	
	Pin to rod coupling	201001-0534	φ5/3-4	
	4-hole pin clamp	201004-0534	φ5/3-4	
	Connecting rod	CF05-250	φ5×250	
	Connecting rod	CF05-200	φ5×200	
	Post	201009-0530	φ5/30°	
	Bone screw	SS20401240	φ4×120	
Humeral shaft fixator GX201009	Compression/Distr action tube	201012-15	φ15	
	Tube to rod coupling	201007-1505	φ15/5	
	4-hole pin clamp	201004-0534	φ5/3-4	
	Post	201009-0530	φ5/30°	
	Bone screw	SS20401240	φ4×120	




Upper limb external fixator instrument sets(A type)














Name	NO.	Size	Picture	
Rod to rod coupling	201002-0505	φ5/5		
Pin to rod coupling	201001-0534	φ5/3-4		
4-hole pin clamp	201004-0534	φ5/3-4		
Peri-articular pin clamp	201003-0534	φ5/3-4		
Post	201008-0500	φ5/0°		
	201009-0530	φ5/30°		
Self drilling/Self tapping bone screws	SS20300830	HA3×80		
	SS20301035	HA3×100		
	SS20400830	HA4×80		
	SS20401035	HA4×100		
	SS20401240	HA4×120		
Guides	QX2007-03	φ3		
	QX2007-04	φ4		
Carbon connecting rod	CF05-120	φ5×120		
	CF05-150	φ5×150		
	CF05-180	φ5×180		
	CF05-200	φ5×200		
T wrench	QX2004-05	φ5×124		
Stabilization/ Reduction wrench	QX2001-05	φ5×125		
Elbow joint movement device	201011-05	φ5×126		
Adaptor	QX2008-03	φ5×127		
	QX2008-04	φ5×128		
Screw driver	QX2005-03	φ5×129		
	QX2005-04	φ5×130		
Hand drill	QX2003	φ5×131		
Thumb wheel	QX2009-57	φ5×132		
Box	QXH01	φ5×133		
Total:				

Combinational external fixator(A type)

















Combination	Name	NO.	Size	Picture
Proximal tibia external fixator GX201101	Rod to rod coupling	201102-0808	φ8/8	
	Pin to rod coupling	201101-0856	φ8/5-6	
	5-hole pin clamp	201103-0856	φ8/5-6	
	Post	201106-0830	φ8/30°	
	Semi-circular curved rod	201107-08180	φ8/180	
	Connecting rod	AL08-250	φ8×250	
Femoral shaft fixator GX201102	Bone screw	SS20501545	φ5×150	
	Rod to rod coupling	201102-0808	φ8/8	
	Pin to rod coupling	201101-0856	φ8/5-6	
	Connecting rod	AL08-200	φ8×200	
Half ring fixator (Standard) GX201103	Connecting rod	AL08-250	φ8×250	
	Bone screw	SS20501545	φ5×150	
	Rod to rod coupling	201102-0808	φ8/8	
	5-hole pin clamp	201103-0856	φ8/5-6	
	Post	201106-0830	φ8/30°	
Proximal femur fixator GX201104	Connecting rod	AL08-300	φ8×300	
	Bone screw	SS20501545	φ5×150	
	Rod to rod coupling	201102-0808	φ8/8	
	5-hole pin clamp	201103-0856	φ8/5-6	
	Post	201106-0830	φ8/30°	
Pelvic fixator GX201105	Connecting rod	AL08-200	φ8×200	
	Bone screw	SS20501545-B	φ5×150	
	Rod to rod coupling	201102-0808	φ8/8	
	5-hole pin clamp	201103-0856	φ8/5-6	
	Post	201106-0830	φ8/30°	
Ankle joint fixator GX201106	Connecting rod	AL08-300	φ8×300	
	Bone screw	SS20501545	φ5×150	
	Rod to rod coupling	201102-0808	φ8/8	
	Pin to rod coupling	201101-0856	φ8/5-6	
Ankle fusion joint fixator GX201108	Connecting rod	AL08-250	φ8×250	
	Bone screw	SS20501545	φ5×150	
	Bone screw	SS20301240	φ3×120	
	Pin to rod coupling	201101-0856	φ8/5-6	
	Pin to rod coupling	201001-0534	φ5/3-4	
Knee joint GX201107	Ankle joint movement device	201113-08	φ8	
	Rod to rod coupling	201102-0808	φ8/8	
	Rod to rod coupling	201101-0856	φ8/5-6	
	Knee joint movement device	201115-08	φ8	
	5-hole pin clamp	201103-0856	φ8/5-6	
	Post	201106-0830	φ8/30°	
	Connecting rod	AL08-300	φ8×300	
Femoral fixator GX201109	Bone screw	SS20501545	φ5×150	
	Tube to rod coupling	201104-2008	φ20/8	
	Compression/distraction tube	201108-20	φ20	
	5-hole pin clamp	201103-0856	φ8/5-6	
GX201110	Post	201106-0830	φ8/30°	
	Bone screw	SS20301240	φ3×120	
GX201110	Micro-movement compression/distraction unilateral external fixator	GX201110	Standard	

Lower limb external fixator instrument sets(A type)




Name	NO.	Size	Picture	
Rod to rod coupling	201102-0808	φ8/8		
Pin to rod coupling	201101-0856	φ8/5-6		
5-hole pin clamp	201103-0856	φ8/5-6		
Post	201105-0800	φ8/0°		
	201106-0830	φ8/30°		
Self drilling/Self tapping bone screws	SS20501240	HA5×120		
	SS20501545	HA5×150		
	SS20501850	HA5×180		
Guides	QX2007-05	φ5		
Carbon connecting rod	CF08-200	φ8×200		
	CF08-250	φ8×250		
	CF08-300	φ8×300		
Semi-circular curved rod	201107-08180	φ8/180		
	201107-08200	φ8/200		
Stabilization/ Reduction wrench	QX2001-07	#7		
Torque wrench	QX2010-07	#7		
Thumb wheel	QX2009-57	#5/7		
Adaptor	QX2008-05	φ5		
Hand drill	QX2003	standard		
Screw driver	QX2005-05	φ5		
T wrench	QX2004-07	#7		
Box	QXH02	for lower limb		

Combinational external fixator(B type)				
Combination	Name	NO.	Size	Picture
Half ring type GX202001	Ringent rod to rod coupling	202003-0808	φ8/8	
	Ringent pin to rod coupling	202004-0845	φ8/4-5	
	Semi-circular curved rod	202012-08	φ8	
	Connecting rod	AL08-300	φ8×300	
	Compression/distract tube	202023-08	φ8	
Condylar fixator GX202002	Bone screw	SS10401240	HA4×120	
	Ringent rod to rod coupling	202003-0808	φ8/8	
	Ringent pin to rod coupling	202004-0845	φ8/4-5	
	Semi-circular curved rod	202012-08	φ8	
	Trapezoidal curved rod	202014-08	φ8	
Tibia joint GX202003	Connecting rod	AL08-250	φ8×250	
	Bone screw	SS10401240	HA4×120	
	Ringent rod to rod coupling	202003-0808	φ8/8	
	Ringent pin to rod coupling	202004-0845	φ8/4-5	
	Semi-circular curved rod	202012-08	φ8	
Z frame (T type) GX202004	Curved rod	202013-08	φ8	
	Connecting rod	AL08-250	φ8×250	
	Bone screw	SS10401240	HA4×120	
	Ringent rod to rod coupling	202003-0808	φ8/8	
	Ringent pin to rod coupling	202004-0845	φ8/4-5	
T frame GX202005	T clamp	202007-0804	φ8/4	
	Connecting rod	CF08-100	φ8×100	
	Connecting rod	CF08-150	φ8×150	
	Bone screw	SS10401240	HA4×120	
	Pin to rod coupling	202010-1156	φ11/5-6	
Straight frame GX202006	T clamp	202007-1105	φ11/5	
	Connecting rod	CF11-300	φ11×300	
	Bone screw	SS10501345	HA5×130	
	Bone screw	SS10501550-B	HB5×150	
	Pin to rod coupling	202010-1156	φ11/5-6	
Z type GX202007	Connecting rod	CF11-200	φ11×200	
	Connecting rod	CF11-150	φ11×150	
	Bone screw	SS10501550-B	HB5×150	
	Pin to rod coupling	202010-1156	φ11/5-6	
	Connecting rod	CF11-200	φ11×200	
T type GX202008	Connecting rod	CF11-150	φ11×150	
	Bone screw	SS10501550-B	HB5×150	
	Pin to rod coupling	202002-1156	φ11/5-6	
	T clamp	202007-1105	φ11/5	
	Connecting rod	CF11-300	φ11×300	
Z frame (Enhancement) GX202009	Bone screw	SS10501345	HA5×130	
	Connecting rod	CF11-200	φ11×200	
	Connecting rod	CF11-150	φ11×150	
	Bone screw	SS10501550-B	HB5×150	
	Rod to rod coupling	202009-1111	φ11/11	
Pelvic fixator GX202011	Pin to rod coupling	202002-1156	φ11/5-6	
	Connecting rod	CF11-200	φ11×200	
	Connecting rod	CF11-150	φ11×150	
	Bone screw	SS10501345	HA5×130	
	Rod to rod coupling	202003-1111	φ11/11	
Straight type GX202010	Ringent rod to rod coupling	202003-1111	φ11/11	
	Ringent pin to rod coupling	202005-1156	φ11/5-6	
	Ringent rod to rod coupling	202003-1108	φ11/8	
	Connecting rod	CF11-200	φ11×200	
	Connecting rod	CF11-300	φ11×300	
Straight type GX202010	Compression/distract tube	202023-08	φ8	
	Bone screw	SS10501550	HA5×150	
	Pin to rod coupling	202002-1156	φ11/5-6	
Straight type GX202010	Connecting rod	CF11-250	φ11×250	
	Bone screw	SS10501345	HA5×130	
	Pin to rod coupling	202002-1156	φ11/5-6	




Parts&instruments of combinational B external fixator

Name	NO.	Size	Picture
Pin to rod(Enhancement)	202001-1105	φ11/5	
Pin to rod(Enhancement)	202002-1156	φ11/5-6	
	202002-1134	φ11/3.5-4	
Ringent rod to rod coupling	202003-1111	φ11/11	
	202003-1108	φ11/8	
	202003-0808	φ8/8	
Ringent pin to rod coupling	202004-0823	φ8/2.5-3.5	
	202004-0845	φ8/4-5	
Ringent pin to rod coupling	202005-1156	φ11/5-6	
	202005-1134	φ11/3-4	
T clamp(Enhancement)	202006-1105	φ11/5	
T clamp	202007-0804	φ8/4	
	202007-1105	φ11/5	
Rod to rod coupling(Enhancement)	202008-1111	φ11/11	
Rod to rod coupling	202009-1111	φ11/11	
	202009-0808	φ8/8	
Pin to rod	202010-0834	φ8/3-4	
	202010-1156	φ11/5-6	
Ringent pin to rod coupling	202011-1111	φ11/11	
Semi-circular curved rod	202012-08	φ8	
Curved rod	202013-08	φ8	
Trapezoidal curved rod	202014-08	φ8	
Pressurizer	QX3008	Standard	
Compression/Distracton rod	202022-08	φ8	
	202023-08	φ8	




Combinational external fixator sets I(B type)

Name	NO.	Size	Picture
Ringent rod to rod coupling	202003-0808	φ8/8	
Ringent pin to rod coupling	202004-0845	φ8/4-5	
Ringent rod to rod coupling	202003-1111	φ11/11	
Ringent pin to rod coupling	202005-1156	φ11/5-6	
T clamp	202007-1105	φ11/5	
Standard screw	SS10401240	HA4×90	
	SS10501135	HA5×110	
	SS10501345	HA5×130	
	SS10601345	HA6×130	
	SS10501550-B	HB5×150	
	SS10601550-B	HB6×150	
Drill bits	Z28200	2.8×200	
	Z38200	3.8×200	
	Z48200	4.8×200	
Screw driver	QX3001-40	φ4	
	QX3001-50	φ5	
Guides	QX3004-40	φ4	
	QX3004-50	φ5	
	QX3004-60	φ6	
Open spanner	QX3003-1010	#10	
Pressurizer	QX3008	Standard	
Connecting rod	CF08-150	φ8×150	
	CF08-200	φ8×200	
	CF08-250	φ8×250	
	CF08-300	φ8×300	
	CF11-150	φ11×150	
	CF11-200	φ11×200	
	CF11-250	φ11×250	
	CF11-300	φ11×300	
	CF11-350	φ11×350	
Boex	QXH03	B- I	

Combinational external fixator sets II(B type)

Name	NO.	Size	Picture	
Pin to rod coupling	202010-1156	φ11/5-6		
Ringent rod to rod coupling	202011-1111	φ11/11		
T clamp	202007-1105	φ11/5		
Bone screws	SS10501135	HA5×110		
	SS10501345	HA5×130		
	SS10501550-B	HB5×150		
	SS10601345	HA6×130		
	SS10601550-B	HB6×150		
Drill bit	Z38200	3.8×200		
	Z48200	4.8×200		
Screw driver	QX3001-50	φ5		
	QX3001-60	φ6		
Guides	QX3004-40	φ4		
	QX3004-50	φ5		
	QX3004-60	φ6		
Open spanner	QX3003-1010	#10		
Pressurizer	QX3008	标准		
Connecting rod	CF11-150	φ11×150		
	CF11-200	φ11×200		
	CF11-250	φ11×250		
	CF11-300	φ11×300		
	CF11-350	φ11×350		
Box	QXH04	B- II		

Combinational external fixator sets III(B type)

Name	NO.	Size	Picture
Pin to rod coupling	2020020834	φ8/3.5-4	
Rod to rod coupling	202009-0808	φ8/8	
Rod to rod coupling	202009-1111	φ11/11	
Pin to rod coupling	202002-1156	φ11/5-6	
T clamp	202007-1105	φ11/5	
Bone screws	SS10400930	HA4×90	
	SS10501135	HA4×110	
	SS10501345	HA5×130	
	SS10601345	HA6×130	
	SS10501550-B	HB5×150	
	SS10601550-B	HB6×150	
Drill bit	Z28200	2.8×200	
	Z38200	3.8×200	
Drill bit	Z48200	4.8×200	
Screw driver	QX3001-40	φ4	
	QX3001-50	φ5	
Guides	QX3004-40	φ4	
	QX3004-50	φ5	
	QX3004-60	φ6	
Pressurizer	QX3008	Standard	
Open spanner	QX3003-1010	#10	
Connecting rod	CF08-150	φ8×150	
	CF08-200	φ8×200	
	CF08-250	φ8×250	
	CF08-300	φ8×300	
	CF11-150	φ11×150	
	CF11-200	φ11×200	
	CF11-250	φ11×250	
	CF11-300	φ11×300	
	CF11-360	φ11×350	
Box	QXH05	B-III	

Get in touch

24/7 professional service

We believe at **EUROTECH** that simplicity is the best way to a successful operation. We design our products to fit the mentioned idea by making life easier for surgeons, nurses and health care professional. We believe that simple product and user-friendly surgical instruments are of the most important factors to get more efficient & reliable outcomes in the theatre.

Beijing Baiou Hengyuan Technology Co., Ltd.

EUROTECH
RELIABLE ORTHOPEDIC
SOLUTION

Aly Thabet
Managing Director

P +86 159 10265446

E aly.thabet@eurotech-surgical.com

EUROTECH

B3215, 3rd Floor, Building 1
No. 18, Jiuxianqiao East Road,
Chaoyang District, Beijing

www.eurotechsurgical.com/

